
# ELECTRONIC SCALE

## WB-100MA/WB-110MA

### INSTRUCTION MANUAL

 Please read this Instruction Manual carefully and keep it handy for future reference.

# Table of Contents

<input type="checkbox"/> Table of Contents	3
<input type="checkbox"/> Applications	4
<input type="checkbox"/> Safety Notes	4
Caution Symbols	4
⚠ Warning	4
⚠ Caution	4
⚠ [Maintenance]	5
⚠ [Instructions for Measurement]	5
<Usage Conditions>	5
<Storage Conditions>	5
<Power Source>	5
<b>1. Product Assembly and Components</b>	<b>6</b>
■ Product Components	6
(Remote Display Version)	
(Column Mounted Version)	
■ Control Panel	7
■ Symbols and their Meanings	7
<b>2. Preparation for Use</b>	<b>8</b>
■ Level Adjustment	8
(Remote Display Version)	
(Column Mounted Version)	
■ Installation of the Equipment	8
(Remote Display Version)	
■ Prepare Power Source	9
<b>3. Operating Instructions</b>	<b>10</b>
1. Standard Weighing Procedures	10
2. Tare Weighing	10
3. Cancellation of Tare Weighing Function	10
4. BMI Calculation:	11
<b>4. RS-232C Communication Specifications</b>	<b>12</b>
■ Specifications	12
■ Names and connection methods for signal lines	12
■ Weighing Mode	12
■ BMI Mode	12
<b>5. Troubleshooting</b>	<b>13</b>
■ Possible Problems	13
<input type="checkbox"/> Specifications	Back cover

## Applications

- This equipment has been calibrated as a precision weighing instrument and can be used to certify weights and/or business transactions.
- It can be used to obtain reference data during medical examinations, such as periodic checkups, and can help the prevention of obesity.

## Safety Notes

### Caution Symbols

In this Instruction Manual and the descriptions on the equipment, there are various caution symbols that will enable you to use this equipment safely and correctly so that harm to you and other people and loss of your property will be prevented.

The symbols and their meanings are as follows. Please understand them well before proceeding with the Instruction Manual.

-  **Warning** This symbol indicates the risk of death or serious injury if the safety instructions are ignored and the equipment is misused.
-  **Caution** This symbol indicates the risk of physical injury and/or damage to the equipment if the safety instructions are ignored and the equipment is misused.
-  This symbol indicates general precautions that should be taken when using this equipment.

### Warning

- To prevent electric shocks,  
Do not insert or remove the power plug with wet hands.
- To prevent fire hazards,  
Always use a correctly wired 230V AC not multiple extension cords.

### Caution

- To prevent the spread of infections,  
Always clean the Weighing platform after use. Sterilize the platform with an appropriate alcohol-based cleaning fluid before each use.
- Installation of the Weighing Platform  
Ensure the weighing platform is placed on a level and stable surface. If the equipment is used when the Platform is unstable, for example because not all the feet are on a level surface, there is a risk that it will topple over and it will be impossible to obtain an accurate measurement.  
To avoid the risk of injury and damage to the equipment, never jump on the platform.
- How to handle the equipment when moving it  
Since the Platform is heavy, please lift it by holding the Platform firmly with both hands so that you do not drop it.
- AC Adapter  
Ensure that the adapter supplied (model A30930G) is always used with this equipment. Other adapters may cause malfunctions.
- Battery  
When changing batteries, ensure that the new batteries are inserted with the correct polarity  $\oplus/\ominus$ . If the polarity is incorrect, then the batteries may leak and damage the equipment.

- The circuit board incorporates a fuse (100mA, 125V) for electronic equipment. If devices other than the original AC adapter are used or a large electrical current is applied, the fuse will be blown for safety purposes. If the fuse is blown, the unit must be inspected and repaired. As the fuse forms a fixed part of the unit and cannot be easily disassembled by the customer, please consult your nearest Tanita sales office or agent when the fuse needs to be replaced.

### [Maintenance]

This is a precision manufactured and accurately calibrated product. Please observe the following instructions.

- Never disassemble or adjust the equipment, as this may cause malfunctions. The equipment should only be inspected in accordance with the regulations in each country.
- When not in use for a long time, unplug the AC Adapter from the wall socket.
- When not in use for a long time, remove the batteries before storing the equipment.
- Do not wipe the equipment with corrosive chemicals (benzine, cleaner, etc.). Please use a neutral detergent to clean the equipment.
- Avoid subjecting the equipment to excessive shocks or vibrations.
- When disposing of this unit, please do so in accordance with the prevailing regulations in each country.
- If an unauthorised person attempts to disassemble or repair any of the part, the warranty will become invalid. When the unit malfunctions, please consult your nearest Tanita sales office or agent.

### [Instructions for Measurement]

Please observe the following instructions for accurate measurement.

- Measurement may not be possible on a surface that is strongly vibrating. In such cases transfer the equipment to a stable surface with little or no vibration.
- When the equipment has been transferred to another location with a temperature difference of more than 20°C, leave at least two hours before using.
- While taking measurements, please keep any persons using transmitters, such as a mobile phone, away from the unit to avoid causing margin errors.

#### <Usage Conditions>

Temperature Range for Use : 0°C / +35°C  
Relative Humidity : 30% / 80% (without condensation)

#### <Storage Conditions>

Temperature Range of Environment : -10°C / +50°C  
Range of Relative Humidity : 10% / 90% (without condensation)  
To avoid malfunctions, do not store the equipment in direct sunlight, significant temperature changes, the risk of dampness, a large amount of dust, in the vicinity of fires or where there is a risk of receiving vibrations or shocks.


#### <Power Source>

	Value Range
Voltage Range	230VAC
Frequency Range	50Hz
Electric Current Range	65mA


# 1.Product Assembly and Components

## Product Components Please verify that the following components are in the package.


(Remote Display Version)

(Column Mounted Version)

## Control Panel

- Power on and/or reset the display to zero. **NET mark:** An arrow ◀ appears next to this icon when the tare weighing function is on.
- Power off. **STABILIZED mark:** Indicates that the displayed figure has stabilized.
- Tare weighing function. **STEP ON mark (Ready mark):** Indicates that the scale is ready to begin measurements.
- BMI mode, and/ or height setting. **HEIGHT mark (Unit of height):** Indicates the height is set
- Mode selection or increase numerical values. **BMI mark:** Indicates that BMI calculation.
- Mode selection or decrease numerical values. **BMI MODE mark:** Indicates when BMI mode is on.
- Weight Lock Release. **WEIGHT LOCK mark:** Indicates when the measurement is stabilized. \* It appears only when the Weight Lock Switch is on.


## Symbols and their Meanings

	: Display Unit On.		: Display Unit Off.		: Direct current
	: Class II Equipment		: Type B applied part		: Zero reset
	: Turn on the tare weighing function.		: Caution. Refer to the attached notes.		: Conformity with Medical Device Directive 93/42/EEC

## 2.Preparation for Use


### Level Adjustment

#### (Remote Display Version)

\*Adjust the level by rotating the adjustable feet.

#### (Column Mounted Version)


\*Adjust the level by rotating the adjustable feet.


For accurate measurement, place the unit on as flat a surface as possible and adjust the level with the adjustable feet so the bubble in the level gauge is in the center of the frame.

### Caution

Place the equipment on a level surface. If the equipment is not stable, for example because not all the feet are on a level surface, there is a risk that it will topple over and it will be impossible to obtain an accurate measurement.


### Installation of the Equipment (Remote Display Version)

#### When using the stand



Please pull out the stand at the back of the display box.  
When putting away the display box, please insert the stand in the slot on the back of the display box.

#### When installing the display box on the wall


- Screw the wall-hanging screws into the wall (leave the screws protruding from the wall by 5 mm).
- Insert the screw heads in the holes on the back of the display box.
- Pull down the display box securely into place.

- If the Display Box is not fixed sufficiently, there is a danger that it may fall down.
- Avoid installing the Display Box in a place where users may step on or get caught by the cable.
- Do not swing the Display Box by the cable.

### Prepare Power Source

#### ① When using the AC Adapter

Insert the cable plug for the AC Adapter into the DC jack on the side of the Display Box, and insert the AC Adapter plug into a 230V AC mains DC jack.

Then confirm the Pilot Lamp on the AC Adapter lights up.


### Warning

- To avoid electric shocks, do not insert and remove the plug with wet hands.

### Caution

- Use only the AC adapter (MODEL: A30930G) specified. Using an AC Adapter other than that specified may cause malfunctions.

#### ② When using dry cell batteries

Unscrew the Fixing Screws for the Battery Compartment on the right side of the Display Box and pull out the Battery Compartment. Place dry cell batteries LR6 (AA) in the correct positions as instructed inside the Battery Compartment and fix the Battery Compartment with the Fixing Screws.

If  $\text{Lo}$  is displayed on the Display, please replace all the six batteries LR6 (AA) with new ones.

$\text{Lo}$  appears when the battery voltage is lower than 6.5 V  $\pm 10\%$ .

- Please change the batteries LR6 (AA) carefully to avoid dropping them on your feet.
- Ensure that the batteries are inserted with the correct polarity  $\oplus/\ominus$ . If the polarity is incorrect, then the batteries may leak and damage the equipment.
- When not in use for a long time, remove the batteries before storing the equipment.
- When the equipment has been transferred to another location with a temperature difference of more than 20°C, leave at least two hours before using.

#### ③ Auto Power off Interval Setting

When the Scale power off, Press  $\text{PIN MODE SET}$  key and hold.

Press  $\text{ON/OFF}$  key.

Release the holding  $\text{PIN MODE SET}$  key.


The display turns to "Auto Power Off Interval Setting Mode" showing "20" as 20 minutes power off setting as default.

Press  $\uparrow$  (up arrow) or  $\downarrow$  (down arrow) key and set 10 for 10 minutes and/or CONT for continuous setting if required.

Press  $\text{PIN MODE SET}$  key to complete.


### 3. Operating Instructions

#### 1. Standard Weighing Procedures


1. Turn on the power by pressing the key.
  2. After all the segments flash, [0.0kg] is displayed.
  3. Step onto the center of the scale platform.  
The weight reading will be displayed.
  4. Measurement is completed.  
When the person being weighed steps off the scale, the display will return to [0.0kg].
- Turn off the power by pressing the key.

#### 2. Tare Weighing


1. Turn on the power by pressing the key.
2. After all the segments flash, [0.0kg] is displayed.
3. Carefully place the tare items (clothes, etc.) on the center of the platform.  
The weight of the items will be displayed.
4. Press the key.  
An arrow will appear next to the mark and [0.0kg] will be displayed.
5. When weighing a person/object including tare items, the weight of the person/object minus that of the tare items, which is the net weight, will be displayed.
  - To change the tare weight, cancel the tare weighing function and repeat from step 3.
  - When using the tare weighing function, the measurement range is reduced by the tare weight.  
(Measurement Range) = (Maximum Measurable Weight) - (Tare Weight)

#### 3. Cancellation of Tare Weighing Function

1. Remove the tare items (clothes, etc.) from the platform.  
The weight will be displayed as a negative value.
2. Press the key.  
The arrow next to the [NET] will disappear and [0.0kg] will be displayed.

#### 4. BMI Calculation:

1. Turn on the power by pressing the key.
  2. After all the segments flash, [0.0kg] is displayed.
  3. Push for within two seconds.  
[0.0kg] and the default height [170cm] will be displayed.
  4. Set your height with the key.  
(This can also be entered after you step on the scale).  
  
Stand still on the scale platform.  
The weight will be displayed.
  5. Once your height and weight have been established, press .  
The weight and BMI figure will be displayed.
- \* By pressing the key in BMI mode, the scale will switch alternately between BMI and height (and can be changed with the key) display.
- \* The display will return to [0.0 kg] and the default height [170cm], if it changes to a height display after stepping off the scale during measuring BMI.

#### Guidelines based on NIH/WHO BMI

Under Weight		BMI < 18.5
Normal Range		18.5 BMI < 25
Pre obese		25 BMI < 30
Obese	CLASS I	30 BMI < 35
Obese	CLASS II	35 BMI < 40
Obese	CLASS III	BMI 40


## 4. RS-232C Communication Specifications

### Specifications

Communication standard	EIA RS-232C
Communication method	Asynchronous method
Signal speed	2,400 baud
Data bit length	7 bit
Parity	EVEN
Stop bit	1 bit

### Names and connection methods for signal lines

Terminal no.	Signal name
1	
2	TXD
3	RXD
4	
5	GND
6	
7	
8	
9	

- ⚠ An RS-232C connector (D sub 9 pin female) is attached to the side of the display box.
- When connecting the display box with a personal computer, etc., please use an RS-232C straight cable.
- Use the RS232C cable shorter than 3m long.

### Weighing Mode

When receiving 4 bytes of data as indicated in 1-1 while using the scale in weighing mode, weight data will be output in 14 bytes, as indicated in 1-2.

#### 1-1 Reception data

1	2	3	4
D	R	CR	LF

#### 1-2 Transmission data

1	2	3	4	5	6	7	8	9	10	11	12	13	14
S	D					7	5	.	0	k	g	CR	LF

Details of Weighing Mode Transmission Data			
Byte nos. 1-2	Header	2 bytes fixed	SD: Unstable S_: Stable SI: Overload, etc.
Byte nos. 3-4	Separator	2 bytes fixed	2 spaces
Byte nos. 5-10	Weight data	6 bytes fixed	Right-aligned
Byte nos. 11-12	Measurement units	2 bytes fixed	Kg:kg
Byte nos. 13-14	Terminator	2 bytes fixed	CR+LF

### BMI Mode

While the scale is used in BMI mode, the output data will vary depending on the form of the received data.

- When receiving 4 bytes of data, as indicated in 1-1, the scale will calculate your BMI using the currently input height figure, and 26 bytes of data will be output, as indicated in 2-1 or 2-2, depending on the measurement unit.
- When receiving 11 bytes of data, including measurement units and height figures, in the form indicated in 1-2 or 1-3, the scale will calculate your BMI using these figures, and 26 bytes of data will be output, as indicated in 2-1 or 2-2, depending on the measurement unit.

#### 1-1 Reception data

1	2	3	4
D	R	CR	LF

#### 1-2 Reception data (height data when using kg)

1	2	3	4	5	6	7	8	9	10	11
D	H	k	g			1	8	0	CR	LF

\*Please input 90-99cm height data as 090-099.

#### 2-1 Transmission data (kg/cm)

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	
S	D					7	5	.	0	k	g	,			1	8	0	,			2	3	.	1	CR	LF

Details of BMI mode transmission data			
Byte nos. 1-2	Header	2 bytes fixed	SD: Unstable S_: Stable SI: Overload, etc.
Byte nos. 3-4	Separator	2 bytes fixed	2 spaces
Byte nos. 5-10	Weight data	6 bytes fixed	Right-aligned
Byte nos. 11-12	Measurement units	2 bytes fixed	Kg:kg
Byte nos. 13	Separator	1 byte fixed	1 comma
Byte nos. 14-18	Height data	5 bytes fixed	Right-aligned
Byte nos. 19	Separator	1 byte fixed	1 comma
Byte nos. 20-24	BMI data	5 bytes fixed	Right-aligned
Byte nos. 25-26	Terminator	2 bytes fixed	CR+LF

## 5. Troubleshooting

### Possible Problems

If problems occur, please refer to the following instructions. If the problem persists contact your nearest Tanita sales or repair center.


Problem	Solution
The unit does not turn on when the  key is pressed.	<ul style="list-style-type: none"> <li>When Using with the AC Adapter Ensure that the AC Adapter cable is properly connected with the Display Unit. Ensure that the AC Adapter is correctly inserted in the mains DC jack.</li> <li>Using dry cell batteries The batteries may have run out. Confirm that the equipment is working by using the AC Adapter or new batteries.</li> </ul>
Measurement is not correct.	<ul style="list-style-type: none"> <li>Please confirm that all the adjustable feet are resting on a level surface.</li> </ul>
Measurement does not become stable. Display turns on and off.	<ul style="list-style-type: none"> <li>Measurement may be impossible in the case of excessive vibration. Use the equipment in another place.</li> </ul>
 is displayed.	<ul style="list-style-type: none"> <li>It is possible that the maximum weight capacity has been exceeded. Please stop measuring.</li> </ul>
 is displayed.	<ul style="list-style-type: none"> <li>Do not stand on or put objects on the Platform before starting measurement. Start measurement after confirming there is nothing on the Platform.</li> </ul>

**TANITA**  
Monitoring Your Health

## Specifications

TYPE		WB-100MA	WB-110MA	
Power Source		AC adapterCentre minusMODEL A30930GCLASS 2 Input Voltage : 230VAC50Hz65mA Output Voltage : 9VDC±0.45VDC Rated Current : 300mA No Load Input Voltage : 12VDC±0.6VDC LR6 (AA Alkaline battery) x 6, (Not included)		
Power Consumption		0.3W max		
Weight Measurement	Measurement System	Strain Gauge Load Cell		
	Maximum Capacity/Minimum Graduation	200kg/0.1kg	270kg/0.1kg	
	Accuracy at first calibration	±0.1kg		
Input Items	Height	90 – 249cm/1cm increments		
Output Items	Display	Weight	200kg/0.1kg increments	270kg/0.1kg increments
		Height	90 – 249cm/1cm increments	
		BMI	0.1 increments	
Size (Remote Display Version)	Display	159×209×56mm		
	Platform	301×336×80mm		
	Cable Length	2m		
Size (Column Mounted Version)		336×600×1070mm		
Display		Upper 5Digits and Lower 4Digits LCD Height of numerals 25mm		
Output Data Interface		RS-232C (D-sub 9pins Female Connector)		
Temperature Range of usage		0°C/+35°C		
Relative Humidity		30%/80% (without condensation)		
Weight of Equipment (Remote Display Version)		5.1kg (except battery)		
Weight of Equipment (Column Mounted Version)		10.6kg (except battery)		
Rated Power		DC9V 300mA (LR6 (AA Alkaline battery) ×6 not included)		
Battery Life		approximately 100 hours of continuous use when using LR6 (AA Alkaline battery)		

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This product meets the following requirements :

1. Medical Device Directive (93/42/EEC)

Safety Standard : EN60601-1:1990,AM No.1:1993,AM No.2:1995  
IEC 60601-1:1988,AM No.1:1991,AM No.2:1995

EMC Standard : EN60601-1-2:2001  
IEC60601-1-2:2001

2.Non-Automatic weighing Instruments (90/384/EEC)

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## TANITA Corporation
## TANITA Corporation of America, inc.
Arlington Heights, Illinois 60005 U.S.A.
## TANITA Health Equipment H.K.LTD.
## TANITA Europe GmbH
## TANITA France S.A.