



*WHO = World Health Organization (OMS: Organisation Mondiale de la Santé)

germanengineering

BM 49 Tensiomètre à Bras - Fonction vocale

Mesure entièrement automatique de la tension artérielle et du pouls au niveau du bras

Avec fonction vocale (4 voix disponibles):

- (1): D, F, I, NL
- (2): GB, ES, PT, GR
- (3): RO, PL, CZ, HU

Volume sonore réglable, voix désactivable

Grand écran très lisible

Moyenne de toutes les mesures enregistrées et de la tension du matin et du soir des 7 derniers jours

2 x 60 emplacements de mémoire

Indicateur OMS* (échelle des valeurs selon les directives de l'OMS*)

Détection de l'arythmie: met en garde en cas de turbulence du rythme cardiaque

Produit médical

Pour un tour de bras de 22 - 36 cm

Date et heure / Arrêt automatique

Message en cas d'erreur d'utilisation

Présélection de pression et déflation

Tension artérielle

Indicateur du niveau des piles

Pochette de rangement pratique

4 piles AAA de 1,5V incluses

Taille (L//H): 57,0 x 47,0 x 25,5 mm

Poids: env. 359 g (Piles incluses)

Pièces disponibles séparément N° d'art.: 071.60

3 ans de garantie

UE: 4 / Carton de transport: 4 x 4

- N° EAN: (1): 4211125656287
 (2): 4211125656294
 (3): 4211125656317

N° d'art.: (1): 656.28 (2): 656.29 (3): 656.31

BM 49 Speaking blood pressure monitor

Fully automatic blood pressure and pulse measurement on the upper arm

With talking function (4 languages adjustable):

- (1): D, F, I, NL
- (2): GB, ES, PT, GR
- (3): RO, PL, CZ, HU

Volume adjustable respectively possibility to turn of speaking function completely

Large, easy-to-read display

Average value stored of all measurements and morning and evening blood pressure for the last 7 days

2 x 60 memory spaces

With WHO* indicator (classification of values according to WHO* guidelines)

Arrhythmia detection: warning in case of possible heart rhythm disturbance

Medical product

For upper arm circumferences from 22 - 36 cm

Date and time /automatic switch-off

Signal in case of application error

Automatic pressure preselection and deflation

Blood pressure passport

Battery changing indicator

Practical storage pouch

Incl. 4 x 1.5V AA batteries

Product measurements (L//W//H): 57.0 x 47.0 x 25.5 mm

Product weight: approx. 359 g (Incl. batteries)

Mains part available separately Item no.: 071.60

3 year guarantee

Sales unit: 4 / Shipping carton: 4 x 4

- EAN-No.: (1): 4211125656287
 (2): 4211125656294
 (3): 4211125656317

Item.-No.: (1): 656.28 (2): 656.29 (3): 656.31